


**DECTRIS**  
detecting the future



# QUADRO


*The fastest direct electron detector for  
Materials Science*





The new QUADRO detector delivers unmatched performances to Materials Science TEM and pushes the limits of electron detection. At last, a direct electron detector capable of single electron counting at very high count-rates with exceptional dynamic range, speed and sensitivity, and no worries to damage the detector even with a focused electron beam.

The QUADRO is powered by our new EIGER2 ASIC, entirely designed by DECTRIS and featuring our patented INSTANT RETRIGGER® technology delivering 24-bits dynamic range and up to 18'000 frames per second readout in windowed mode, no binning, no gimmicks!



*Diffraction Pattern from an Evaporated Aluminum standard calibration sample by TedPella.*

*The dynamic range of QUADRO allows to resolve the feeble rings at higher angle (above  $0.82\text{\AA}^{-1}$ ) without the need of masking the unscattered beam.*



*Collect all the data your sample delivers!*

### Key Advantages

- Direct electron detection
- Region of interest feature
- Up to 18'000 frames/sec
- No dead time
- Noise-free electron counting
- $10^7$  counts/pixel/sec
- No beam-stop required
- Ideal DQE also at low energies

### Applications

- Electron diffraction
- 4D-STEM
- Strain mapping
- Lorentz microscopy mapping
- Ptychography
- *In-Situ* TEM
- Dynamic TEM
- LEEM/PEEM




High resolution imaging and diffraction of a single monolayer of Molybdenum Disulfide.

The exagonal lattice is clearly resolved in imaging mode, even with a low-end microscope, and its diffraction pattern required only  $\sim 2 \cdot 10^5$  electrons and 10 ms of acquisition time.

QUADRO is faster than the instabilities in your lab. The high speed and sensitivity allow the acquisition of high-quality data in spite of drifts and interferences.

Upgrade your microscope, not your room!




## Technical specifications

### QUADRO

Pixel size [ $\mu\text{m}^2$ ]	75 x 75
Number of pixels	514 x 514
Active area, width x height [ $\text{mm}^2$ ]	38.6 x 38.6
Energy range [kV]	30 - 300
Threshold range [kV]	10 - 80
Maximum frame rate, ROI [Hz]	18'000
Maximum frame rate, full frame [Hz]	2'250 (16-bit), 4'500 (8-bit)
Sensor material	Si or CdTe
Data format	HDF5
Cooling	Water, 20 °C
Detector dimensions (WHD) [ $\text{mm}^3$ ]	200 x 350 x 200
Detective Quantum Efficiency (0)	0.99 @ 100 kV; 0.96 @ 200 kV
Detector mounting	Bottom, on-axis
Radiation hardness [el./ $\text{mm}^2$ ]	$>5 \times 10^{15}$

All specifications are subject to change without notice.



**DECTRIS Ltd.**  
Täferweg 1  
5405 Baden-Dättwil  
Switzerland

Phone +41 56 500 21 00  
Fax +41 56 500 21 01  
E-Mail [sales@dectris.com](mailto:sales@dectris.com)  
[www.dectris.com](http://www.dectris.com)

Registered trademarks: „DECTRIS“: EU, JP, CN, KR (IR0911969), USA (5,253,168)  
“detecting the future“: EU, CN, KR, JP, AUS (1191333), USA (4,607,800)  
“DECTRIS Instant Retrigger“: EU, CN, KR, JP, AUS (1224728), USA (4,797,363)  
“DECTRIS QUADRO“: CH (722 238), CN (IR1437238)  
© 2018 DECTRIS Ltd., all rights reserved • Subject to technical modifications;  
Rev.0. • Printed 05/2019