


**DECTRIS**  
detecting the future



## DECTRIS ELA

A new hybrid pixel  
electron detector



### Technical specifications

Frame rate	2,250 Hz at 16 bits
	4,500 Hz at 8 bits
	18,000 Hz with ROI mode
Count rate capability	> 1 pA/pixel
Number of pixels	1,024 x 512
Sensor material	Silicon
Image bit-depth	32 bits
Energy range	30-200 keV
Point-spread function	< 1.3 pixel at 200 keV


All specifications are subject to change without notice.

### New era of EELS with DECTRIS ELA


Powered by DECTRIS Hybrid Pixel Technology, ELA is the ideal electron-counting detector for Electron Energy Loss Spectroscopy (EELS), offering unprecedented speed, sensitivity, and dynamic range. It can handle probe currents well over 100 pA, simultaneously recording the intense zero-loss peak and the core levels as well as counting single electrons in all regions of the spectrum. Its high frame rate allows for fast elemental mapping in one go, which is especially crucial when working with beam-sensitive materials. Thanks to its dead-time free readout, you will not lose electrons between images, collecting all the information that your sample can deliver. Combined with a high-resolution spectrometer, the DECTRIS ELA detector will significantly enhance applications such as EELS and 4D-STEM, allowing your experiments to be conducted in minutes. It will increase instrument throughput and reduces operation time significantly.



Plotkin-Swing, Benjamin, et al. „Hybrid pixel direct detector for electron energy loss spectroscopy.“ *Ultramicroscopy* (2020): 113067.



Fast elemental mapping of an STO/BTO/LMSO multilayer, including the unsaturated ZLP, acquired as 32 separate 128 x 128 spectrum images of 8 s each, aligned and summed (4.3 min total). From Left to Right: ZLP intensity map, similar to a bright-field image. Map of Ti L<sub>2,3</sub> edge. Map of Ba M<sub>4,5</sub> edge. Map of La M<sub>4,5</sub> edge. RGB composite map. (Image courtesy: NION Ltd.)



EELS of hexagonal boron nitride with  $E_0 = 60$  keV, beam current of 105 pA, 100 s exposure, and energy dispersion set to 0.5 eV/channel, acquired using a Nion Iris EEL spectrometer on a Nion UltraSTEM 200MC monochromated scanning transmission electron microscope. (Image courtesy: NION Ltd.)