

EIGER Geometry

Throughout the present document we use the notations from DectrisGeometryDocumentation. This document defines the rotation matrix \mathbf{R} and the translation vector \mathbf{t} , which transform the position \mathbf{P}' in the detector system to its equivalent position \mathbf{P} in the lab system through:

$$P_i = R_{ij}P'_j + t_i \quad (1)$$

In addition we introduce the coordinates of the beam centre in the detector system $\mathbf{c} = (c_0, c_1, 0)$ and the distance between the detector and the sample d . \mathbf{c} and d are related to \mathbf{R} and \mathbf{t} through

$$R_{ij}c_j + t_i = \delta_{i3}d \quad (2)$$

where $\delta_{ij} = 1$ if $i = j$ else 0.

Given \mathbf{R} and \mathbf{t} (2) defines \mathbf{c} and d , which is not true vice versa. So the EIGER API adopts the convention that whenever \mathbf{c} or d changes, \mathbf{R} is kept constant and \mathbf{t} is computed from (2):

$$t_0 = -R_{00}c_0 - R_{01}c_1 \quad (3)$$

$$t_1 = -R_{10}c_0 - R_{11}c_1 \quad (4)$$

$$t_2 = R_{20}c_0 + R_{21}c_1 - d \quad (5)$$

Whenever \mathbf{R} or \mathbf{t} changes \mathbf{c} and d is computed from (2):

$$c_0 = \frac{R_{01}t_1 - R_{11}t_0}{R_{00}R_{11} - R_{01}R_{10}} \quad (6)$$

$$c_1 = \frac{R_{10}t_0 - R_{00}t_1}{R_{00}R_{11} - R_{01}R_{10}} \quad (7)$$

$$d = R_{20}c_0 + R_{21}c_1 - t_2 \quad (8)$$