


Technical Specification and Operating Procedure

PILATUS 100K-S Detector System



Version 1.8



Table of Contents

1	DOCUMENT HISTORY	3
1.1	CHANGES.....	3
2	HOW TO USE THIS TECHNICAL SPECIFICATION	4
2.1	ADDRESS AND SUPPORT	4
2.2	EXPLANATION OF SYMBOLS	5
2.3	EXPLANATION OF TERMS	5
2.4	USE OF THE PILATUS 100K-S	5
3	TECHNICAL SPECIFICATION	6
3.1	RATINGS.....	7
3.2	AMBIENT CONDITIONS	8
4	DIMENSIONS AND CONNECTORS	9
4.1	PILATUS 100K-S.....	10
4.1.1	<i>Front side of the Detector</i>	10
4.1.2	<i>Backside of the Detector</i>	10
4.1.3	<i>The status LEDs</i>	11
4.1.4	<i>Connectors and Connecting Cables/Pipes</i>	11
4.2	POWER SUPPLY.....	13
4.3	COMPUTER	14
5	INSTALLING THE DETECTOR SYSTEM	15
5.1	MOUNTING.....	15
5.1.1	<i>Mounting from Above</i>	15
5.1.2	<i>Mounting from the Bottom</i>	16
5.2	GROUNDING OF THE DETECTOR SYSTEM	17
5.3	CONNECTION TO NITROGEN OR DRY AIR.....	18
5.4	CONNECTING THE CABLES	19
6	TEMPERATURE AND HUMIDITY CONTROL	20
7	GETTING STARTED	21
7.1	STARTUP SEQUENCE	21
7.2	FIRST COMMANDS.....	22
8	TURNING OFF THE DETECTOR	23
9	STORING THE DETECTOR.....	23
10	CLEANING AND MAINTENANCE	24
11	FAULTS.....	25
12	CERTIFICATION TESTS	26



1 Document History

Actual document

Version	Date	status	prepared	checked	released
1.8	29.08.2014	released	PH	SB	SB

1.1 Changes

Version	Date	Changes	released
1.0	05.11.2008	First release	PS
1.1	17.04.2009	Various improvements	PS
1.4	29.01.2010	Connection cables and temperature sensors	BS
1.4.1	18.02.2010	Signal levels	BS
1.5	09.04.2010	Dimensions and grounding	BS
1.6	04.10.2010	Pictures	BS
1.7	22.07.2011	Thicker sensors and conformity with standards	BS
1.8	29.07.2014	Adapted Detector External Trigger Input specification	SB



2 How to use this Technical Specification

Before you start to operate the PILATUS 100K-S detector system please read this technical specification and the user manual thoroughly.

The technical specification and the user manual together form the user documentation.

2.1 Address and Support

DECTRIS Ltd.
Neuenhoferstrasse 107
5400 Baden
Switzerland
Phone: +41 56 500 21 00
Fax: + 41 56 500 21 01

Email: support@dectris.com






Should you have questions concerning the system or its use, please contact us via phone, mail or fax.



Before you ship the system back, please contact us to receive the necessary transport and shipping information.



2.2 Explanation of Symbols

Symbol	Description
	Important or helpful notice
	Caution. Please follow the instructions carefully to prevent equipment damage or personal injury.
	DC-current
	AC-current
	Ground

2.3 Explanation of Terms

Term	Description
MCB	Module Control Board
DCB	Detector Control Board
DAC	Digital to Analog Converter

2.4 Use of the PILATUS 100K-S

The PILATUS 100K-S detector system has been designed for the detection of X-rays from synchrotrons or laboratory sources.

It is intended for indoor use only.

For other applications, please contact DECTRIS for additional information.



Do not use the detector in vacuum

The PC can be mounted in a standard 19 inch rack, which has to be properly grounded.



Make sure that the PC has adequate ventilation.




3 Technical Specification

Number of modules	1
Sensor	Reverse-biased silicon diode array
Sensor thickness	320 μm 450 μm
Quantum efficiency	3 keV: 48% 3 keV: 48% 8 keV: 95% 8 keV: 96% 15 keV: 51% 15 keV: 64%
Pixel size	172 x 172 μm^2
Module size	83.8 x 33.5 mm^2
Format	487 x 195 = 94'965 pixels
Dynamic range	20 Bits = 1'048'576
Counting rate per pixel	$> 2 \times 10^6$ X-ray/sec
Energy range	4.5 – 36 keV
Energy resolution	500 eV
Adjustable threshold range	4 – 18 keV
Threshold dispersion	50 eV
Readout time	2.3 ms
Framing rate	300 Hz
Point-spread function	1 pixel
Data formats	Raw data, TIF, EDF, CBF
External trigger/gate	3.3 - 5V, 3 different modes
Software interface	Through socket connection; Clients for EPICS, SPEC and stand-alone operation are available
Cooling	Air-cooled
Operating temperature (internal)	20-35 °C
Dimensions (W x H x D)	143 x 110 x 180 mm^3
Weight	Approx. 2.3 kg



3.1 Ratings

Device	Definition
Detector Power Input	+12 V DC, 1.5 A, 18 Watt
Power supply	100 – 240 VAC, 50/60 Hz
Detector External Trigger Input	2.1 V - 5.0 V High level 0.0 V - 0.8 V Low level 5.0 V absolute maximum  Applying a higher voltage will destroy the input.
Detector Enable output	5 V TTL (max. 100 mA)
PC	100 - 240 VAC, 6.5 A 50/60 Hz 345 Watt Can be connected to all common supply voltages.



3.2 Ambient Conditions

The PILATUS 100K-S detector is designed only for indoor use. The following ambient conditions must be fulfilled:

Condition	Range
Operating temperature	+20° to +35°C
Operating humidity	< 80% at 20°C, non-condensing
Storage temperature	+15° to +40°C
Storage humidity	< 40% at 20°C, non-condensing



Operate the detector only with connected nitrogen or dry air supply. For details see section 5.3.



When storing the detector make sure the temperature and humidity inside the transport box doesn't exceed the specified range. Use drying agent.



If the detector system is stored at low temperature, make sure that no condensation moisture develops.



The PILATUS 100K-S detector is equipped with a temperature and humidity control, see section 6.



4 Dimensions and Connectors




Figure 1. Drawing of the PILATUS 100K-S detector (printed separately in the user documentation folder)



4.1 Pilatus 100K-S

4.1.1 Front side of the Detector

The detector comes with a protective cover (1 mm, low carbon steel sheet metal 1.0330, St12) for the front window which should be removed for operation. The sensors are behind a 12 μm thick Mylar[®] (PET) foil coated with 100 nm aluminum to protect it from dust and touch.



Do not touch the Mylar[®] foil.



The cover may not protect the detector from direct synchrotron beam.



Figure 2. PILATUS 100K detector with cover in place (left) and removed (right).

4.1.2 Backside of the Detector




Figure 3. PILATUS 100K detector viewed from the back. With cable cover mounted (left) and dismantled (right).





4.1.3 The status LEDs

LED	Description
Power	If green, all supply voltages are ok. If red, module power is off. This can be caused by over/under temperature, too much humidity or one of the internal voltages failed.
EN	Yellow, detector is making an exposure.

4.1.4 Connectors and Connecting Cables/Pipes

Connector	Description
DATA	<p>Data connection</p> <p>Cable: Use only the included cable. Custom made cables or connectors will not work.</p> <p> Minimum bending radius in a fix installation: 5 times the cable diameter. Minimum bending radius in a flexible installation: 15 times the cable diameter.</p> <p>Details of pins:</p> <div style="display: flex; align-items: center; justify-content: center;"> <div style="border: 1px solid black; padding: 5px; margin-right: 20px;"> <p>Data in/out 9 Pol Sub-D</p> </div> <div style="display: flex; flex-direction: column; align-items: center;"> <div style="width: 10px; height: 20px; background-color: black; margin-bottom: 10px;"></div> <div style="width: 10px; height: 20px; background-color: black;"></div> </div> <div style="display: flex; flex-direction: column; align-items: center;"> <div style="width: 10px; height: 20px; background-color: black; margin-bottom: 10px;"></div> <div style="width: 10px; height: 20px; background-color: black;"></div> </div> </div> <div style="margin-top: 10px;"> <p>Pin 1: Rx In Pin 5: Tx Out Pin 6: Rx In (GND) Pin 9: Tx Out (GND)</p> </div> <div style="margin-top: 10px;"> <div style="border: 1px solid black; padding: 5px; margin-right: 20px;"> <p>PCI-Karte Rx 9 Pol Sub-D Pin 5 in, Pin 9 GND</p> </div> <div style="border: 1px solid black; padding: 5px;"> <p>PCI-Karte Tx 9 Pol Sub-D Pin 1 in, Pin 6 GND</p> </div> </div>
Power	Main voltage 12 VDC from external power supply
EXT IN	<p>External Trigger Input</p> <p>Use a Lemo ® Type 00 (NIM-CAMAC) cable.</p> <p>For ratings see 3.1.</p>
EN OUT	<p>TTL output signal; high when counting is enabled.</p> <p>Use a Lemo ® Type 00 (NIM-CAMAC) cable.</p> <p>For ratings see 3.1.</p>



Connector	Description
A OUT	Analog output signal; to monitor the analog signal of a pixel. Use a Lemo ® Type 00 (NIM-CAMAC) cable.
Ground 	Functional ground of the detector system.  Although the detector may be grounded via the mounting bolts, the detector can be grounded via its ground connector to establish a defined grounding.
N2	Nitrogen for humidity control. For details see section 5.3. Pipe: Use a pipe with outer diameter of 4 mm.

The fixation pieces are stored at the bottom of the detector.



4.2 Power Supply

The power supply is a compact switching power unit.



Figure 4. Power supply +12V DC 60W

Connector	Description
Mains Input	100V -240VAC 50/60Hz 700mA Can be connected to all common supply voltages.
Output	+12V DC / 2.5A



4.3 Computer

The computer is a high power server with a proprietary data acquisition card to communicate with the detector.



The PC should be kept behind a firewall and should not have outside internet access.



The operating system is optimized for high speed data acquisition and has a custom kernel: Therefore, do not permit any software upgrades on the kernel!



Do not install or run any other software on the computer, except as tools and software necessary to configure your data acquisition protocol.



Figure 5. Connectors at the PCI Mezzanine Card for data acquisition

Connector	Description
RX	Receive data
TX	Transmit data



5 Installing the Detector System

5.1 Mounting

The detector can be mounted in any position and has 4 holes with M6 threads in the base plate. There are two ways to mount the detector:

5.1.1 Mounting from Above

Use the detachable mounting pieces, which are stored on the base plate of the detector.



Make sure the mounting pieces are mounted and properly tightened




Figure 6. Mounting the detector, bottom view

Mounting points when mounting from the top



5.1.2 Mounting from the Bottom

From the bottom with four M6 bolts directly into the base plate




Figure 7. Mounting the detector, bottom view

Mounting points when mounting from the bottom



10 mm.

The four M6 bolts should not intrude into the detector more than



Make sure the detector is properly mounted.



5.2 Grounding of the Detector System



The main plug of the computer has to be connected to a grounded power outlet.



Although the detector might be grounded via the mounting bolts, the detector can be grounded additionally via the functional ground connector at the back (M4 screw-in tap hole) to establish a defined grounding.



Figure 8. Ground connector on the back side of the detector



5.3 Connection to Nitrogen or Dry Air

The PILATUS detector has to be connected to a nitrogen or dry air flow to avoid humidity and condensation when it is outside the storage box. For storage see section 9.



Figure 9. Nitrogen connector on the back of the detector



Humidity can damage the detector. Make sure that the detector is operated with connected nitrogen or dry air.

Recommended flow of nitrogen	0.5-1 liter/min.
As alternative	Oil free, dry air of < 2% relative humidity can be used. Recommended flow: 0.5-1 liter/min.
Gas pressure	Minimum 1 bar Maximum 2 bar



5.4 Connecting the Cables

To operate the detector, the data cables and the ground should be connected. For specification of cables and pipes see section 4.1.4.

The PILATUS 100K-S detector is equipped with a data cable.

Connect RX to RX and TX to TX on the GigaSTaR Card in the PC. Labels for the connections are printed on the server housing. The single connectors go to the corresponding connector on the back of the detector.



The data cable should be pulled onto the computer connectors with the screws, rather than forcefully pushed on.

A forceful connection can damage the PCI card.



It is important for data integrity that the screws are tightened.



To plug or unplug any cables, turn the detector off.



6 Temperature and Humidity Control

The PILATUS 100K-S detector head has 2 combined temperature and humidity sensors.

The temperature and humidity control shuts down the power of the detector modules when the humidity or the temperature of the different sensors exceeds or falls below the following limits:

Channel	Location	Shutdown Temperature °C		Shutdown Humidity %
		Low limit	High limit	
0	Power board	15	55	80%
1	Base plate	15	35	80%

The fan and the communication with the PC will remain active after a temperature shut down (only shuts down the power of the modules).



Connect the PILATUS 100K detector to nitrogen or dry air flow. Keep a constant flow rate (according to section 5.3).



7 Getting Started

Before operating the detector, make sure you have read the previous chapters in the technical specification and the user manual.



Check these items before turning the detector system on:

- Mount the detector properly
- Connect the detector to power
- Connect the detector to nitrogen or dry air at the recommended flow rate
- Connect the PC and the detector data cables
- Attach monitor and keyboard to the computer

7.1 Startup Sequence

- Turn on nitrogen or dry air flow if possible at least 30 minutes before turning on the detector
- Turn on the PC
- Start a shell
- The default path is: /home/det
- Change the directory to: p2_det
- Turn on the power switch at the back-side of the detector
- Type runtvx.
Runtvx starts a script file which initializes the detector system and opens the Camserver and TVX windows.



7.2 First commands

See the detailed description of all commands in the User Manual.

Type the following commands in TVX:

- *rbd*; self test of the detector (digital part of the pixel)
- *calibdet*; self test of the detector (analog part of the pixel)
- *cam setthreshold 5900* (example for a *setthreshold* of 5.9 keV): sets the energy threshold of the detector. It is important that a threshold is set. Otherwise the detector is not trimmed.
In normal operation, the threshold should always be set to 50% of the energy of the incoming X-rays. See user manual for more information!
- *expose 10*; creates an image with an exposure time of 10 seconds



8 Turning off the Detector

- Turn off the power switch of the detector
- Do not remove nitrogen/dry air connection and leave it at the recommended flow rate according to section 5.3.

9 Storing the Detector

Even if the detector is not in operation, it is recommended to maintain the nitrogen or dry air flow.

In case the detector is stored, please follow these instructions:

For storage of up to 1 week:

Store the detector in the storage box and add 200 g drying agent (i.e. silica gel) into the storage box.

For storage longer than 1 week:

Pack the detector into a plastic bag, add drying agent into the plastic bag and seal the plastic bag.

Then place the detector in the storage box and add 200 g drying agent (i.e. silica gel) to the storage box.

Check the humidity inside the box frequently to comply with the storage requirements in chapter 3.2



10 Cleaning and maintenance

The housing can be cleaned with a soft tissue.



The Mylar ® foil should not be touched or cleaned

The PILATUS 100K-S detector system is maintenance free.



11 Faults

Problems	Causes	Solution
PC doesn't start properly	<ul style="list-style-type: none">• PC is not properly powered up• PCI card is not properly mounted	<ul style="list-style-type: none">• Depending on the type of the PC, there are switches on the back and on the front of the PC which have to be in the correct position.• Open the PC and check that the PCI cards are properly mounted and tightened
Detector shuts down	Temperature or humidity error See section 6 for shut down values. Check the LEDs on the back of the detector: If the temperature and the power LED is red, a temperature or humidity error has occurred	<ul style="list-style-type: none">• Wait until the detector cools down• Restart the detector again.• Check the temperature of the detector with the command in camserver: type "thread"• Check that the ventilator of the detector is running properly and the ventilation holes are not covered• Check the flow of nitrogen or dry air
Detector shuts down	Over current Threshold level set too low and the detector starts oscillating	<ul style="list-style-type: none">• Increase the threshold level: In TVX type: setthreshold midg 5000
Images look strange after initialization	Detector is not properly initialized	Run the following commands in TVX: <ul style="list-style-type: none">• <i>setdac</i>• <i>calibdet</i>• <i>expose 1</i>
After a trigger command has been issued camserver reports an error	After issuing the trigger command in camserver a trigger has to occur within 15 seconds	Make sure the trigger command occurs within 15 seconds
Detector housing is humid	Ambient humidity around the detector exceeds the operating conditions	Shut down the detector immediately and check the humidity. Power up the detector only when the ambient humidity has been reduced



12 Certification Tests

The product is in conformity with the listed standards.

Applied standards / Normes appliquées / Verwendete Normen

EN 61000-6-2: 2005 (IEC 61000-6-2: 2005)	Electromagnetic compatibility (EMC) - Part 6-2: Generic standards - Immunity for industrial environments Compatibilité électromagnétique (CEM) - Partie 6-2: Normes génériques - Immunité pour les environnements industriels Elektromagnetische Verträglichkeit (EMV) - Teil 6-2: Fachgrundnormen - Störfestigkeit - Industriebereich
EN 61000-6-4: 2007 (IEC 61000-6-4: 2006)	Electromagnetic compatibility (EMC) - Part 6-4: Generic standards - Emission standard for industrial environments Compatibilité électromagnétique (CEM) - Partie 6-4: Normes génériques - Norme sur l'émission pour les environnements industriels Elektromagnetische Verträglichkeit (EMV) - Teil 6-4: Fachgrundnormen - Fachgrundnorm Störaussendung - Industriebereich
EN 61326-1: 2006 (IEC 61326-1: 2005)	Electrical equipment for measurement, control and laboratory use - EMC requirements -- Part 1: General requirements Matériel électrique de mesure, de commande et de laboratoire - Exigences relatives à la CEM -- Partie 1: Exigences générales Elektrische Mess-, Steuer-, Regel- und Laborgeräte - EMV-Anforderungen -- Teil 1: Allgemeine Anforderungen